Bijlage 17 Cardiogen informatie

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CardioGen-82®

Rubidium Rb 82 Generator

HIGHLIGHTS OF PRESCRIBING INFORMATION
 These highlights do not include all the information needed to use **CARDIOGEN-82** safely and effectively. See full prescribing information for **CARDIOGEN-82** at www.fda.gov/medwatch.

CARDIOGEN-82 (rubidium Rb 82 generator)
 To produce rubidium Rb 82 chloride injection, for intravenous use
 Initial U.S. Approval: 1989

WARNING: HIGH LEVEL RADIATION EXPOSURE WITH USE OF INCORRECT ELUENT AND FAILURE TO FOLLOW THE ELUATE TESTING PROTOCOL
 Please see full prescribing information for complete boxed warning

High Level Radiation Exposure with Use of Incorrect Eluent
 Using the incorrect eluent can cause high Strontium (Sr) 82 and Sr 85 breakthrough levels (5-1).

- Use only additive-free 0.9% Sodium Chloride Injection USP to elute the generator (2.5)
- Immediately stop the patient infusion and permanently discontinue the use of the affected CardioGen-82 generator if the incorrect solution is used to elute the generator (4)
- Evaluate the patient's radiation absorbed dose and monitor for the effects of radiation to critical organs such as bone marrow (2.10)

Excess Radiation Exposure with Failure to Follow the Eluate Testing Protocol
 Excess radiation exposure occurs when the levels of Sr 82 or Sr 85 in the rubidium Rb 82 chloride injection exceed limits (5.2)

- Record eluate volume, including waste and test volumes (2.5)
- Strictly adhere to the generator eluate testing protocol (2.6, 2.7)
- Stop using the generator if it reaches any of its Expiration Limits (2.8)

RECENT MAJOR CHANGES----- 10/2020
 Directions for Eluting Rubidium Rb 82 Chloride Injection (2.5)
 CardioGen-82 Infusion System Model 1700 Eluate Testing Protocol (2.7)
 CardioGen-82 Dose Delivery Limit (2.9)

INDICATIONS AND USAGE-----
 CardioGen-82 is a closed system used to produce rubidium Rb 82 chloride injection for intravenous use. Rubidium Rb 82 chloride injection is a radioactive diagnostic agent indicated for Positron Emission Tomography (PET) imaging of the myocardial under rest or pharmacologic stress conditions to evaluate regional myocardial perfusion in adult patients with suspected or existing coronary artery disease. (1)

DOSEAGE AND ADMINISTRATION-----
 • Use only additive free 0.9% Sodium Chloride Injection USP to elute the generator (2.5)
 • Use CardioGen-82 only with the CardioGen-82 Infusion System Model 510 or Model 1700. (2.4)

FULL PRESCRIBING INFORMATION: CONTENTS
WARNING: HIGH LEVEL RADIATION EXPOSURE WITH USE OF INCORRECT ELUENT AND FAILURE TO FOLLOW THE ELUATE TESTING PROTOCOL
1 INDICATIONS AND USAGE
2 DOSAGE AND ADMINISTRATION
 2.1 Radiation Safety – Drug Handling
 2.2 Rubidium Rb 82 Chloride Injection Dosage
 2.3 Image Acquisition Guidelines
 2.4 Infusion System
 2.5 Directions for Eluting Rubidium Rb 82 Chloride Injection
 2.6 CardioGen-82 Infusion System Model 510 Eluate Testing Protocol
 2.7 CardioGen-82 Infusion System Model 1700 Eluate Testing Protocol
 2.8 CardioGen-82 Expiration
 2.9 CardioGen-82 Dose Delivery Limit
3 CONTRAINDICATIONS
4 DOSAGE FORMS AND STRENGTHS
5 WARNINGS AND PRECAUTIONS
 5.1 High Level Radiation Exposure with Use of Incorrect Eluent
 5.2 Excess Radiation Exposure with Failure to Follow the Eluate Testing Protocol
 5.3 Risks Associated with Pharmacologic Stress
 5.4 Volume Overload
 5.5 Cumulative Radiation Exposure: Long-Term Risk of Cancer
6 ADVERSE REACTIONS
 6.1 Postmarketing Experience

FULL PRESCRIBING INFORMATION
WARNING: HIGH LEVEL RADIATION EXPOSURE WITH USE OF INCORRECT ELUENT AND FAILURE TO FOLLOW THE ELUATE TESTING PROTOCOL
High Level Radiation Exposure with Use of Incorrect Eluent
 Patients are exposed to high radiation levels when the CardioGen-82 generator is eluted with the incorrect eluent due to high Sr 82 and Sr 85 breakthrough levels. See Warnings and Precautions (5.1).
 • Use only additive-free 0.9% Sodium Chloride Injection USP to elute the generator. (See Dosage and Administration (2.5))
 • Immediately stop the patient infusion and permanently discontinue the use of the affected CardioGen-82 generator if the incorrect solution is used to elute the generator. (See Contraindications (4)).
 • Evaluate the patient's radiation absorbed dose and monitor for the effects of radiation to critical organs such as bone marrow. (See Dosage and Administration (2.10))
Excess Radiation Exposure with Failure to Follow the Eluate Testing Protocol
 Excess radiation exposure occurs when the levels of Sr 82 or Sr 85 in the rubidium Rb 82 chloride injection exceed specified limits. (See Warnings and Precautions (5.2))
 • Record each generator eluate volume, including waste and test volumes, and keep a record of the cumulative eluate volume. (See Dosage and Administration (2.5))
 • Strictly adhere to the generator eluate testing protocol, to minimize the risk of excess radiation exposure, including daily testing and additional testing at Alert Limits. (See Dosage and Administration (2.6, 2.7))
 • Stop using the generator if it reaches any of its Expiration Limits.
 o An eluate Sr 82 level of 0.01 µCi/mL Rb 82, or
 o An eluate Sr 85 level of 0.1 µCi/mL Rb 82. (See Dosage and Administration (2.8))

INDICATIONS AND USAGE

CardioGen-82 is a closed system used to produce rubidium Rb 82 chloride injection for intravenous administration. Rubidium Rb 82 chloride injection is indicated for Positron Emission Tomography (PET) imaging of the myocardium under rest or pharmacologic stress conditions to evaluate regional myocardial perfusion in adult patients with suspected or existing coronary artery disease.

DOSEAGE AND ADMINISTRATION

2.1 Radiation Safety - Drug Handling

- Rubidium Rb 82 is a radioactive drug and should be handled with appropriate safety measures to minimize radiation exposure during administration. (See Warnings and Precautions (5.5))
- Use only additive-free 0.9% Sodium Chloride Injection USP to elute the generator. (See Boxed Warning, Contraindications (4), and Warnings and Precautions (5.1))
- Limit the use of radiopharmaceuticals to physicians who are qualified by training and experience in the safe use and handling of radionuclides and whose experience and training have been approved by the appropriate government agency authorized to license the use of radionuclides.
- Wear waterproof gloves and effective shielding when handling rubidium Rb 82 chloride injection and the infusion system.
- Observe aseptic techniques in all drug handling, including rubidium Rb 82 chloride injection and the infusion system.
- Visually inspect the drug for particulate matter and discoloration prior to administration whenever solution and container permit. Do not administer eluate from the generator if there is any evidence of foreign matter.

2.2 Rubidium Rb 82 Chloride Injection Dosage

The recommended adult single dose of rubidium Rb 82 chloride injection per rest or stress component of a PET myocardial perfusion imaging (MPI) procedure is 1,480 MBq (40 mCi) with a range of 1,110 MBq to 2,220 MBq (30 mCi to 60 mCi).

- Do not exceed a single dose of 2,220 MBq (60 mCi).
- Use the lowest dose necessary, consistent with the goal of as low as reasonably achievable (ALARA), to obtain adequate cardiac visualization.
- The choice of dosage and infusion rate should be made at the discretion of the physician performing the Rb 82 PET myocardial perfusion imaging scan, based on his/her consideration of multiple factors, including patient body mass, and the imaging equipment and acquisition methodology used to perform the procedure. For example, 3D image acquisition may require doses at the lower end of the recommended range, compared to 2D imaging.
- Administer the single dose at 50 mL/minute, or at 20 mL/minute (CardioGen-82 Infusion System Model 1700 only), through a catheter inserted into a large peripheral vein; do not exceed a total infusion volume of 100 mL.
- Instruct patients to void as soon as a study is completed and as often as possible thereafter for at least one hour.
- The maximum available activity (delivery limit) will decrease as the generator ages. (See Dosage and Administration (2.9)).

2.3 Image Acquisition Guidelines

Administer two separate single doses to complete rest and stress myocardial perfusion imaging as follows:
 For Rest Imaging:
 • Administer a single ("rest") rubidium Rb 82 chloride dose;
 • Start imaging 60 to 90 seconds after completion of the infusion of the rest dose, and acquire images for 5 minutes.
 For Stress Imaging:
 • Administer a pharmacologic stress agent in accordance with its prescribing information;
 • Administer the pharmacologic stress agent, administer the second dose of Rb 82 at the time interval according to the prescribing information of the pharmacologic stress agent;
 • Start imaging 60 to 90 seconds after completion of the stress Rb 82 chloride dose infusion, and acquire images for 5 minutes.

For Both Rest and Stress Imaging:
 • If a longer circulation time is anticipated (e.g., in a patient with severe left ventricular dysfunction), start imaging 120 seconds after the rest dose.
 • Acquisition may be started immediately post-injection if dynamic imaging is needed.

2.4 Infusion System

Use CardioGen-82 only with the CardioGen-82 Infusion System Model 510 or the CardioGen-82 Infusion System Model 1700, which are specifically designed for use with the generator and capable of accurate measurement and delivery of doses of rubidium Rb 82 chloride injection.

- If using the CardioGen-82 Infusion System Model 510, refer to the Eluate Testing Protocol in Section 2.6. (See Dosage and Administration (2.7)).
- Follow instructions in the CardioGen-82 Infusion System Model 510 or Model 1700 Operator's Manual for the set up and intravenous infusion of rubidium Rb 82 chloride injection doses(1,2)
- The generator, when used with either infusion system, provides ± 10% accuracy for rubidium Rb 82 chloride injection between doses of 1,110 MBq to 2,220 MBq (30 mCi to 60 mCi).

2.5 Directions for Eluting Rubidium Rb 82 Chloride Injection

- Use only additive-free 0.9% Sodium Chloride Injection USP to elute the generator. (See Boxed Warning, Contraindications (4), and Warnings and Precautions (5.1)).
- Prepare the 0.9% Sodium Chloride Injection USP for use with the Saline Bag.
- Prepare the intravenous port in accordance with the DOSAGE AND ADMINISTRATION section of the approved prescribing information of the 0.9% Sodium Chloride Injection USP.
- The intravenous administration port of the sodium chloride container must be penetrated only one time.
- Strap the saline bag provided with the CardioGen-82 infusion system on the additive-free 0.9% Sodium Chloride Injection USP container and install on the CardioGen-82 Infusion System.
- Once the container port closure is penetrated, it should remain installed on the CardioGen-82 Infusion System for its entire period of use.
- Allow at least 10 minutes between elutions for regeneration of Rb 82.
- If the CardioGen-82 Infusion System Model 510 is used, discard the first 50 mL eluate each day the generator is eluted, and employ proper safety precautions since the eluate contains radioactivity. If the CardioGen-82 Infusion System Model 1700 is used, the system will automatically discard the first 50 mL of eluate each day the generator is eluted.
- If the CardioGen-82 Infusion System Model 510 is used, maintain an on-going record of all eluate volumes (washing, testing, dosing volumes) including a summary of the cumulative volume of eluate from the generator. If the CardioGen-82 Infusion System Model 1700 is used, the system software automatically records and saves all eluate volumes (all flushing, Quality Control (QC) testing, and patient infusions), representing the cumulative volume of eluate each day the generator is eluted.

2.6 CardioGen-82 Infusion System Model 510 Eluate Testing Protocol

Use only additive-free 0.9% Sodium Chloride Injection USP for all elutions. (See Boxed Warning, Contraindications (4), and Warnings and Precautions (5.1)).

Observe aseptic technique throughout.

Follow all instructions in the CardioGen-82 Infusion System Model 510 Operator's Manual for performing all eluate testing as described. Before administering rubidium Rb 82 chloride injection to the first patient each day, perform the following testing:

- Strontium Alert Limits and Mandatory Eluate Testing:
 o Use an ionization chamber-type dose calibrator for eluate testing.
 o Daily, before administering rubidium Rb 82 chloride injection to any patient, perform an eluate testing to determine Rb 82, Sr 82, and Sr 85 levels.
 o Perform additional daily eluate tests after detecting any of the following Alert Limits:
 o 14 L total eluate volume has passed through the generator column, or
 o Sr 85 level reaches 0.002 µCi per mCi Rb 82, or
 o Sr 85 level reaches 0.02 µCi per mCi Rb 82.
 o For example, if an Alert Limit were reached and the clinical site eluted less than 750 mL from the generator during the day, then no additional eluate tests would be performed that day.
- o If the same clinical site the next day eluted 1,500 mL from the generator, then the site would have performed three tests that day: 1) the required daily test that precedes any patient dosing, 2) a test at the 750 mL elution point, and 3) a test at the 1,500 mL elution points.
- o If a generator's Alert Limit is reached, the clinical site performs the additional daily tests (at intervals of 750 mL) each subsequent day the generator is used. The additional tests are necessary to promptly detect excessive Sr 82 and/or Sr 85 in eluates.

Rubidium Eluate Level Testing:

1. Set a dose calibrator for Rb 82 as recommended by the manufacturer or use the Co-60 setting and divide the reading obtained by 0.548. Obtain the reading from the instrument.
2. Elute the generator with 50 mL of Saline. Chloride Injection USP and discard the eluate (first elution).
3. Allow at least 10 minutes for the eluate to stoppage in the syringe. Note the exact time of end of elution (E.O.E.).
4. 50 mL into and collect the eluate in a stoppered glass vial (plastic containers are not suitable). Note the exact time of end of elution (E.O.E.).
5. Using the dose calibrator, determine the activity of Rb 82 and note the time of the reading. Correct the reading for decay to the E.O.E. using the appropriate decay factor for Rb 82 (see Table 1). Note: If the reading is taken 2 ½ minutes after end of elution, multiply the dose calibrator reading by 4 to correct for decay.

Strontium Eluate Level Testing:
 5. Using the sample obtained for the Rb 82 activity determination, allow the sample to stand for at least one hour to allow for the complete decay of Rb 82.
 6. Measure the activity of the sample in a dose calibrator at the setting recommended by the manufacturer for Rb 82 and/or Sr 82. As an alternative, use the Co-60 setting and the reading obtained divided by 0.548. **Set the instrument to read in microcuries and record in the display.**
 7. Calculate the ratio (R) of Sr 85/Sr 82 on the day (postcalibration) of the measurement using the ratio of Sr 85/Sr 82 on the day of calibration provided on the generator label and the Sr 85/Sr 82 Ratio Factor from Table 2. Determine R using the following equation:

$$R = \frac{[Sr\ 85]}{[Sr\ 82]} \times \text{calibration date X Ratio Factor on the day (post calibration) of measurement}$$

8. Use a correction factor (F) of 0.478 to compensate for the contribution of Sr 85 to the reading.
9. Calculate the amount of Sr 82 in the sample using the following equation:

$$Sr\ 82\ (\mu Ci) = \frac{[Sr\ 82]}{\text{dose calibration reading } (\mu Ci)} \times (1 + (R) \cdot F)$$

Example: dose calibrator reading (µCi) = 0.8; Sr85/Sr82 ratio (R) = 1.48; correction factor (F) = 0.478.

$$Sr\ 82\ (\mu Ci) = \frac{0.8}{\text{dose calibration reading } (\mu Ci)} \times (1 + (1.48) \cdot 0.478)$$

10. Determine if Sr 82 in the eluate exceeds an Alert or Expiration Limit by multiplying the result obtained in step 10 by (R) of Rb 82 at End of Elution (see below for further instructions based on the Sr 82 level)
 Example: 0.47 µCi or Sr 82: 50 mCi of Rb 82 E.O.E.
 0.47 µCi Sr 82 _____ = 0.0094 µCi/mCi Rb 82 (is above Alert Limit of 0.002; additional daily eluate testing must be performed)

11. Determine if Sr 85 in the eluate exceeds an Alert or Expiration Limit by multiplying the result obtained in step 10 by (R) as calculated in step 7 (above).
 Example: 0.0094 x 1.48 = 0.014 µCi Sr 85/mCi Rb 82 (test result is below Alert and Expiration Limits)

Use Table 1 to calculate the decay factor for Rb 82; step 4 (above).

TABLE 1		
Physical Decay Chart: Rb 82 half-life 75 seconds		
Seconds	Fraction Remaining	Fraction Remaining
0*	1.000	0.218
15	0.871	0.190
30	0.758	0.165
45	0.660	0.144
60	0.574	0.125
75	0.500	0.109
90	0.435	0.095
105	0.379	0.083
120	0.330	0.072
135	0.287	0.063
150	0.250	

*Elution time

Use Table 2 to calculate the ratio (R) of Sr 85/Sr 82; step 7 (above).

TABLE 2		
Sr 85/Sr 82 Ratio Chart (Sr 85 T½ = 65 days, Sr 82 T½ = 25 days)		
Days	Ratio Factor	Days
0*	1.00	15
1	1.02	16
2	1.03	17
3	1.05	18
4	1.07	19
5	1.09	20
6	1.11	21
7	1.13	22
8	1.15	23
9	1.17	24
10	1.19	25
11	1.21	26
12	1.23	27
13	1.25	28
14	1.27	29

*Day of calibration

2.7 CardioGen-82 Infusion System Model 1700 Eluate Testing Protocol

Use only additive-free 0.9% Sodium Chloride Injection USP for all elutions. (See Boxed Warning, Contraindications (4), and Warnings and Precautions (5.1)).

Observe aseptic technique throughout.

Follow all instructions in the CardioGen-82 Infusion System Model 1700 Operator's Manual for performing all eluate testing as described. Before administering rubidium Rb 82 chloride injection to the first patient each day, perform the following testing:

- Strontium Alert Limits and Mandatory Eluate Testing:
 o Use an ionization chamber-type dose calibrator for eluate testing.
 o Daily, before administering rubidium Rb 82 chloride injection to any patient, perform an eluate testing to determine Rb 82, Sr 82, and Sr 85 levels.
 o Perform additional daily eluate tests after detecting any of the following Alert Limits:
 o 14 L total eluate volume has passed through the generator column, or
 o Sr 85 level reaches 0.002 µCi per mCi Rb 82, or
 o Sr 85 level reaches 0.02 µCi per mCi Rb 82.
 o For example, if an Alert Limit were reached and the clinical site eluted less than 750 mL from the generator during the day, then no additional eluate tests would be performed that day.
- o If the same clinical site the next day eluted 1,500 mL from the generator, then the site would have performed three tests that day: 1) the required daily test that precedes any patient dosing, 2) a test at the 750 mL elution point, and 3) a test at the 1,500 mL elution points.
- o If a generator's Alert Limit is reached, the clinical site performs the additional daily tests (at intervals of 750 mL) each subsequent day the generator is used. The additional tests are necessary to promptly detect excessive Sr 82 and/or Sr 85 in eluates.

Rubidium Eluate Level Testing:

1. Set a dose calibrator for Rb 82 as recommended by the manufacturer or use the Co-60 setting and divide the reading obtained by 0.548. Obtain the reading from the instrument.
2. Elute the generator with 50 mL of Saline. Chloride Injection USP and discard the eluate (first elution).
3. Allow at least 10 minutes for the eluate to stoppage in the syringe. Note the exact time of end of elution (E.O.E.).
4. 50 mL into and collect the eluate in a stoppered glass vial (plastic containers are not suitable). Note the exact time of end of elution (E.O.E.).
5. Using the dose calibrator, determine the activity of Rb 82 and note the time of the reading. Correct the reading for decay to the E.O.E. using the appropriate decay factor for Rb 82 (see Table 1). Note: If the reading is taken 2 ½ minutes after end of elution, multiply the dose calibrator reading by 4 to correct for decay.

Strontium Eluate Level Testing:
 5. Using the sample obtained for the Rb 82 activity determination, allow the sample to stand for at least one hour to allow for the complete decay of Rb 82.
 6. Measure the activity of the sample in a dose calibrator at the setting recommended by the manufacturer for Rb 82 and/or Sr 82. As an alternative, use the Co-60 setting and the reading obtained divided by 0.548. **Set the instrument to read in microcuries and record in the display.**
 7. Calculate the ratio (R) of Sr 85/Sr 82 on the day (postcalibration) of the measurement using the ratio of Sr 85/Sr 82 on the day of calibration provided on the generator label and the Sr 85/Sr 82 Ratio Factor from Table 2. Determine R using the following equation:

$$R = \frac{[Sr\ 85]}{[Sr\ 82]} \times \text{calibration date X Ratio Factor on the day (post calibration) of measurement}$$

8. Use a correction factor (F) of 0.478 to compensate for the contribution of Sr 85 to the reading.
9. Calculate the amount of Sr 82 in the sample using the following equation:

$$Sr\ 82\ (\mu Ci) = \frac{[Sr\ 82]}{\text{dose calibration reading } (\mu Ci)} \times (1 + (R) \cdot F)$$

Example: dose calibrator reading (µCi) = 0.8; Sr85/Sr82 ratio (R) = 1.48; correction factor (F) = 0.478.

$$Sr\ 82\ (\mu Ci) = \frac{0.8}{\text{dose calibration reading } (\mu Ci)} \times (1 + (1.48) \cdot 0.478)$$

$$Sr\ 82\ (\mu Ci) = \frac{0.8}{[1 + (1.48)(0.478)]} = 0.47$$

10. Determine if Sr 82 in the eluate exceeds an Alert or Expiration Limit by multiplying the result obtained in step 10 by (R) of Rb 82 at End of Elution (see below for further instructions based on the Sr 82 level)

Example: 0.47 µCi or Sr 82: 50 mCi of Rb 82 E.O.E.
 0.47 µCi Sr 82 _____ = 0.0094 µCi/mCi Rb 82 (is above Alert Limit of 0.002; additional daily eluate testing must be performed)

11. Determine if Sr 85 in the eluate exceeds an Alert or Expiration Limit by multiplying the result obtained in step 10 by (R) as calculated in step 7 (above).
 Example: 0.0094 x 1.48 = 0.014 µCi Sr 85/mCi Rb 82 (test result is below Alert and Expiration Limits)

Use Table 1 to calculate the decay factor for Rb 82; step 4 (above).

TABLE 1		
Physical Decay Chart: Rb 82 half-life 75 seconds		
Seconds	Fraction Remaining	Fraction Remaining
0*	1.000	0.218
15	0.871	0.190
30	0.758	0.165
45	0.660	0.144
60	0.574	0.125
75	0.500	0.109
90	0.435	0.095
105	0.379	0.083
120	0.330	0.072
135	0.287	0.063
150	0.250	

*Elution time

Use Table 2 to calculate the ratio (R) of Sr 85/Sr 82; step 7 (above).

TABLE 2		
Sr 85/Sr 82 Ratio Chart (Sr 85 T½ = 65 days, Sr 82 T½ = 25 days)		
Days	Ratio Factor	Days
0*	1.00	15
1	1.02	16
2	1.03	17
3	1.05	18
4	1.07	19
5	1.09	20
6	1.11	21
7	1.13	22
8	1.15	23
9	1.17	24
10	1.19	25
11	1.21	26
12	1.23	27
13	1.25	28
14	1.27	29

*Day of calibration

2.7 CardioGen-82 Infusion System Model 1700 Eluate Testing Protocol

Use only additive-free 0.9% Sodium Chloride Injection USP for all elutions. (See Boxed Warning, Contraindications (4), and Warnings and Precautions (5.1)).

Observe aseptic technique throughout.

Follow all instructions in the CardioGen-82 Infusion System Model 1700 Operator's Manual for performing all eluate testing as described. Before administering rubidium Rb 82 chloride injection to the first patient each day, perform the following testing:

- Strontium Alert Limits and Mandatory Eluate Testing:
 o Use an ionization chamber-type dose calibrator for eluate testing.
 o Daily, before administering rubidium Rb 82 chloride injection to any patient, perform an eluate testing to determine Rb 82, Sr 82, and Sr 85 levels.
 o Perform additional daily eluate tests after detecting any of the following Alert Limits:
 o 14 L total eluate volume has passed through the generator column, or
 o Sr 85 level reaches 0.002 µCi per mCi Rb 82, or
 o Sr 85 level reaches 0.02 µCi per mCi Rb 82.
 o For example, if an Alert Limit were reached and the clinical site eluted less than 750 mL from the generator during the day, then no additional eluate tests would be performed that day.
- o If the same clinical site the next day eluted 1,500 mL from the generator, then the site would have performed three tests that day: 1) the required daily test that precedes any patient dosing, 2) a test at the 750 mL elution point, and 3) a test at the 1,500 mL elution points.
- o If a generator's Alert Limit is reached, the clinical site performs the additional daily tests (at intervals of 750 mL) each subsequent day the generator is used. The additional tests are necessary to promptly detect excessive Sr 82 and/or Sr 85 in eluates.

Rubidium Eluate Level Testing:

1. Set a dose calibrator for Rb 82 as recommended by the manufacturer or use the Co-60 setting and divide the reading obtained by 0.548. Obtain the reading from the instrument.
2. Elute the generator with 50 mL of Saline. Chloride Injection USP and discard the eluate (first elution).
3. Allow at least 10 minutes for the eluate to stoppage in the syringe. Note the exact time of end of elution (E.O.E.).
4. 50 mL into and collect the eluate in a stoppered glass vial (plastic containers are not suitable). Note the exact time of end of elution (E.O.E.).
5. Using the dose calibrator, determine the activity of Rb 82 and note the time of the reading. Correct the reading for decay to the E.O.E. using the appropriate decay factor for Rb 82 (see Table 1). Note: If the reading is taken 2 ½ minutes after end of

Immediately stop the patient infusion and discontinue use of the affected CardioGen-82 generator if the incorrect eluent is used and evaluate the patient's radiation absorbed dose and monitor for the effects of radiation to critical organs such as bone marrow. When solutions containing calcium ions are used to elute the generator, high levels of radioactivity are present in any subsequent eluate, even with the use of additive free 0.9% Sodium Chloride Injection USP. *[See Boxed Warning, Dosage and Administration (2.10) and Contraindications (4)].*

5.2 Excess Radiation Exposure with Failure to Follow Eluate Testing Protocol
 Excess radiation exposure occurs when the Sr 82 and Sr 85 levels in rubidium Rb 82 chloride injections exceed the specified generator eluate limits. Strictly adhere to the eluate testing protocol to minimize radiation exposure to the patient. Stop using the rubidium generator when the expiration limits are reached. *[See Dosage and Administration (2.6, 2.7) and (2.8)].*

5.3 Risk Associated with Pharmacologic Stress
 Pharmacologic induction of cardiovascular stress may be associated with serious adverse reactions such as myocardial infarction, arrhythmia, hypotension, bronchoconstriction, and cerebrovascular events. Perform pharmacologic stress testing in accordance with the pharmacologic stress agents prescribing information and only in the setting where cardiac resuscitation equipment and trained staff are readily available.

5.4 Volume Overload
 Patients with congestive heart failure or the elderly may experience a transitory increase in circulatory volume load.

5.5 Cumulative Radiation Exposure - Long-Term Risk of Cancer
 Rubidium Rb 82 chloride injection, similar to other radiopharmaceuticals, contributes to a patient's overall long-term cumulative radiation exposure. Long-term cumulative radiation exposure is associated with an increased risk of cancer. Use the lowest dose of rubidium Rb 82 chloride injection necessary for imaging and ensure safe handling to protect the patient and health care worker. *[See Dosage and Administration (2.1) and (2.2)].* Encourage patients to void as soon as a study is completed and as often as possible thereafter for at least one hour.

6 ADVERSE REACTIONS
6.1 Postmarketing Experience
 The following serious adverse reactions have been identified during postapproval use of CardioGen-82. Because these reactions are reported voluntarily from a population of uncertain size, it is not always possible to reliably estimate their frequency or establish a causal relationship to drug exposure.

Radiation Exposure
 High level radiation exposure to the bone marrow has occurred in some patients due to Sr 82 and Sr 85 breakthrough in the eluate when an incorrect solution was used to elute the rubidium Rb 82 generator. *[See Boxed Warning and Warnings (5.1)].*
 Excess radiation exposure has occurred in some patients who received rubidium Rb 82 chloride injections at clinical sites where generator eluate testing appeared insufficient. *[See Boxed Warning, Warnings and Precautions (5.2), and Dosage and Administration (2.6 or 2.7)].*

8 USE IN SPECIFIC POPULATIONS

8.1 Pregnancy
 Risk Summary
 There are no data available on the use of rubidium Rb 82 chloride in pregnant women. Animal reproductive studies have not been conducted with rubidium Rb 82 chloride injection. However, all radiopharmaceuticals have the potential to cause fetal harm depending on the fetal stage of development and the magnitude of the radiation dose. If considering rubidium Rb 82 chloride injection administration to a pregnant woman, inform the patient about the potential for adverse pregnancy outcomes based on the radiation dose from rubidium Rb 82 and the gestational timing of exposure. The estimated background risk of major birth defects and miscarriage for the indicated population is unknown. In the U.S. general population, the estimated background risk of major birth defects and miscarriage in clinically recognized pregnancies is 2-4% and 15-20%, respectively.

8.2 Lactation
 Risk Summary
 There is no information regarding the presence of Rb 82 chloride in human milk, the effects on the breastfed infant or the effects on milk production. Due to the short half-life of rubidium Rb 82 (75 seconds), exposure of a breastfed infant through breast milk can be minimized by temporary discontinuation of breastfeeding. *[See Clinical Considerations].* The developmental and health benefits of breastfeeding should be considered along with the mother's clinical need for Rb 82 chloride injection, any potential adverse effects on the breastfed child from Rb 82 or from the underlying maternal condition.

8.3 Minimizing Exposure
 Clinical Considerations
 Exposure to Rb 82 chloride through breast milk can be minimized if breastfeeding is discontinued when Rb 82 chloride injection is administered. Do not resume breastfeeding until at least one hour after completion of rubidium Rb 82 chloride injection infusion.

8.4 Pediatric Use
 Rubidium Rb 82 chloride injection safety and effectiveness in pediatric patients have not been established.

8.5 Geriatric Use
 In elderly patients with a clinically important decrease in cardiac function, lengthen the delay between infusion and image acquisition. *[See Dosage and Administration (2.3)].* Observe for the possibility of fluid overload. *[See Warnings and Precautions (5.4)].*

8.6 Renal Impairment
 Reductions in renal function are not anticipated to alter clearance of rubidium Rb 82 chloride injection because Rb 82 decays to stable Kr-82 with a half-life of 75 seconds and Kr-82 is exhaled through the lungs.

8.7 Hepatic Impairment
 Reductions in hepatic function are not anticipated to alter clearance of rubidium Rb 82 chloride injection.

11 DESCRIPTION
 CardioGen-82 contains accelerator-produced Sr 82 adsorbed on stannic oxide in a lead-shielded column and provides a means for obtaining sterile nonpyrogenic solutions of rubidium Rb 82 chloride injection. The chemical form of Rb 82 is ⁸²RbClO₄. The amount (millicuries) of Rb 82 obtained in each elution will depend on the potency of the generator. When eluted at a rate of 50 mL/minute, each generator eluate at the end of elution should not contain more than 0.02 microcurie of Sr 82 and not more than 0.2 microcurie of Sr 85 per millicurie of rubidium Rb 82 chloride injection, and not more than 1 microgram of tin per mL of eluate.

11.2 Physical Characteristics
 Rb 82 decays by positron emission and associated gamma emission with a physical half-life of 75 seconds.⁴ Table 5 shows the annihilation photons released following positron emission which are useful for detection and imaging studies.

The decay modes of Rb 82 are: 95.5% by positron emission, resulting in the production of annihilation radiation, i.e., two 511 keV gamma rays; and 4.5% by electron capture, resulting in the emission of "prompt" gamma rays of predominantly 776.5 keV. Both decay modes lead directly to the formation of stable Kr-82.⁴

TABLE 5 Principal Radiation Emission Data

Radiation	Mean Percent Per Disintegration	Mean Energy (keV)
Annihilation photons (2)	191.01	511 (each)
Gamma rays	13-15	776.5

The specific gamma ray constant for Rb 82 is 6.1 R/hour-millicurie at 1 centimeter. The first half-value layer is 0.7 centimeter of lead (Pb). Table 6 shows a range of values for the relative attenuation of the radiation emitted by this radionuclide that results from interposition of various thicknesses of lead.⁵ For example, the use of a 7.0 centimeter thickness of Pb will attenuate the radiation emitted by a factor of about 1,000.

TABLE 6 Radiation Attenuation by Lead Shielding

Shield Thickness (Pb) cm	Attenuation Factor
0.7	0.5
2.3	10 ¹
4.7	10 ²
7.0	10 ³
9.3	10 ⁴

4. Strontium Level Test and Dose Constancy:
 a. The infusion system will elute the generator with 50 mL of additive-free 0.9% Sodium Chloride Injection USP at a rate of 50 mL/min into the stoppered vial (which is specifically provided for use with the CardioGen-82 Infusion System Model 1700).
 b. The Sr detector measures the Rb 82 and strontium in the 50 mL elution.
 c. The infusion system software will automatically calculate the Sr 82 and Sr 85 levels on the day (post calibration) of the measurement using the ratio of Sr 85/Sr 82 on the day of calibration provided on the generator label, and using the full exponential decay calculation for each, accounting for the generator's age.
 d. Using the Rb 82 and strontium measurements, the infusion system software will automatically calculate µCi Sr 82/mCi Rb 82 and µCi Sr 85/mCi Rb 82. The GUI will automatically indicate if the results exceed Alert or Expiration Limits.

5. Constancy Check of the Sr detector: The infusion system GUI will prompt the user to perform the constancy check of the Sr detector.
 a. Place the external constancy source in the detector well of the infusion system.
 b. The infusion system software will automatically calculate the constancy of the Sr detector versus the external constancy source when instructed.

2.8 CardioGen-82 Expiration
 If using the CardioGen-82 Infusion System Model 510, stop use of the CardioGen-82 generator once any one of the following Expiration Limits is reached. If using the CardioGen-82 Infusion System Model 1700, the software will automatically indicate, and will stop use of, the CardioGen-82 generator once any one of the following Expiration Limits is reached.
 • A total elution volume of 17 L has passed through the generator column, or
 • 42 days post calibration date, or
 • An eluate Sr 82 level of 0.01 µCi/mCi Rb 82, or
 • An eluate Sr 85 level of 0.1 µCi/mCi Rb 82.

2.9 CardioGen-82 Dose Delivery Limit
 The maximum available activity (delivery limit) will decrease as the generator ages. Certain doses, including the maximum recommended dose (60 mCi (2,220 MBq)) are not achievable for the entire shelf-life of the generator. Table 3 provides an estimate of the maximum available activity of Rubidium Rb 82 (Delivery Limit) as a function of generator age.

Table 3 Rubidium Rb 82 Dose Delivery Limit Based on Generator Age¹

Generator Age (days) ²	Maximum Rubidium Dose (Delivery Limit)
0-17	60 mCi (2,220 MBq)
24	50 mCi (1,850 MBq)
32	40 mCi (1,480 MBq)
42	30 mCi (1,110 MBq)

¹ Estimate is based on a 100 mCi (3,700 MBq) Sr 82 generator at calibration.
² Generator age at which delivery limit is reached varies with generator activity at release. For example, a 90 mCi (3,330 MBq) generator and a 150 mCi (5,550 MBq) generator will reach a delivery limit <60 mCi at ≥ 14 days and ≥ 33 days, respectively.

2.10 Radiation Dosimetry
 The estimated absorbed radiation doses for Rb 82, Sr 82, and Sr 85 from an intravenous injection of rubidium Rb 82 chloride are shown in Table 4.

Table 4 Adult Absorbed Radiation Dose Coefficient

Organ ^{a,b}	Rb 82 (Average for Rest and Stress) (µSv/3.7MBq) ^c	Sr 82 (µSv/3.7MBq) ^c	Sr 85 (mrem/µCi (µSv/3.7MBq)) ^c
Adrenals	7.56	10.6	5.03
Bone - Osteogenic cells Bone Surface	1.86	107	9.81
Brain	0.60	8.29	2.96
Breast	0.82	7.03	1.72
Gall Bladder/Wall	3.17	8.47	2.82
Heart Wall	16.5	8.18	2.67
Kidneys	20.04	9.18	2.50
Liver	4.20	8.10	2.50
Lower Large Intestine Wall	2.84	5.18	5.14
Lungs ^d	10.7	8.25	2.84
Muscles	1.29	8.14	2.66
Ovaries	1.41	10.2	4.29
Pancreas	8.85	9.10	3.46
Red Marrow	1.19	9.10	9.84
Skin	1.14	7.03	1.75
Small Intestine	4.76	9.62	4.03
Spleen	6.61	8.10	2.54
Stomach	8.14	7.84	2.26
Testes	0.82	7.25	1.70
Thyroid	1.49	8.07	2.33
Thyroid	6.11	6.11	2.57
Upper Large Intestine	5.94	23.7	3.62
Urinary Bladder Wall	1.61	2.19	2.90
Uterus	3.72	9.14	3.32
Total Body	1.77	Not calculated	Not calculated
Effective Dose ^e	4.74 ^f	23.4	4.03

^a Rb 82 doses are averages of rest and stress dosimetry data (see Senthimzichehian et al.^{1,2}). To calculate organ doses (mrem) from Rb 82, multiply the dose coefficient for each organ by the administered activity in mCi.
^b Sr 82 and Sr 85 doses are calculated using software package DCAI and ICRP dose coefficients. To calculate organ doses (mrem) attributable to Sr 82, and Sr 85, multiply the dose coefficients by the calculated amounts of strontium in µCi.³
^c To convert to SI units, insert the dose coefficient into the formula in parentheses, e.g. for adrenals 7.56 mrem/mCi = 7.56 µSv/37 MBq = 2.04 × 10⁻³ Sv/Bq.
^d Calculated from ICRP 66
^e Calculated from ICRP 60
^f Stress phase only

3 DOSAGE FORMS AND STRENGTHS
 CardioGen-82 is a closed system used to produce rubidium Rb 82 chloride injection for intravenous use. CardioGen-82 consists of strontium Sr 82 adsorbed on a hydrous stannic oxide column with an activity of 3,330 MBq (90 millicuries) Sr 82 at calibration time.

4 CONTRAINDICATIONS
 CardioGen-82 is contraindicated if a solution other than additive free 0.9% Sodium Chloride Injection USP has been used to elute the generator at any time. Immediately stop the patient infusion and permanently discontinue the use of the affected CardioGen-82 generator whenever the incorrect eluent is used. *[See Boxed Warning, Contraindications (4), and Warnings and Precautions (5.1)].*

5 WARNINGS AND PRECAUTIONS

5.1 High Level Radiation Exposure with Use of Incorrect Eluent
 Use only additive free 0.9% Sodium Chloride Injection USP to elute the generator. Apply the provided saline tag to the additive free 0.9% Sodium Chloride Injection USP container before use. Additives present in other solutions (particularly calcium ions) expose patients to high levels of radiation by causing the release of large amounts of Sr 82 and Sr 85 into the eluate regardless of the generator's age or prior use. *[See Dosage and Administration (2.1, 2.5, and 2.6, 2.7)].*

Sr 82 (half-life of 25 days (600 hrs.)) decays to Rb 82. To correct for physical decay of Sr 82, Table 7 shows the fractions that remain at selected intervals after the time of calibration.

TABLE 7 Physical Decay Chart: Sr 82 half-life 25 days

Days	Fraction Remaining	Days	Fraction Remaining
0*	1.000	15	0.660
1	0.973	16	0.642
2	0.946	17	0.624
3	0.920	18	0.607
4	0.895	19	0.591
5	0.871	20	0.574
6	0.847	21	0.559
7	0.824	22	0.543
8	0.801	23	0.529
9	0.779	24	0.514
10	0.758	25	0.500
11	0.737	26	0.486
12	0.717	27	0.473
13	0.697	28	0.460
14	0.678	29	0.448

*Calibration time
 To correct for physical decay of Rb 82, Table 1 shows the fraction of Rb 82 remaining in all 15 second intervals up to 300 seconds after time of calibration. *[See Dosage and Administration (2.6, 2.7)].*

12 CLINICAL PHARMACOLOGY

12.1 Mechanism of Action
 Rb 82 is analogous to potassium ion (K⁺) in its biochemical behavior and is rapidly extracted by the myocardium proportional to the blood flow. Rb⁺ participates in the sodium-potassium (Na⁺/K⁺) ion exchange pumps that are present in cell membranes. The intracellular uptake of Rb 82 requires maintenance of ionic gradient across cell membranes. Rb 82 radioactivity is increased in viable myocardium reflecting intracellular retention, while the tracer is cleared rapidly from necrotic or infarcted tissue.

12.2 Pharmacokinetics
 In human studies, myocardial activity was noted within the first minute after peripheral intravenous injection of Rb 82. When areas of infarction or ischemia are present in the myocardium, they are visualized within 2-7 minutes after injection as photon-deficient, or "cold", areas on the myocardial scan. In patients with reduced cardiac function, transit of the injected dose from the peripheral infusion site to the myocardium may be delayed. *[See Dosage and Administration (2.3)].*

Blood flow brings Rb 82 to all areas of the body during the first pass of circulation. Accordingly, visible uptake is also observed in other highly vascularized organs, such as the kidneys, liver, spleen and lungs.

12.3 Pharmacokinetics
 With a physical half-life of 75 seconds, Rb 82 is very rapidly converted by radioactive decay into a trace amount of stable Kr-82 gas, which is passively expired by the lungs.

13 NONCLINICAL TOXICOLOGY

13.1 Carcinogenesis, Mutagenesis, Impairment of Fertility
 No long-term studies have been performed to evaluate carcinogenic potential, mutagenicity potential, or to determine whether rubidium Rb 82 chloride injection may affect fertility in males or females.

14 CLINICAL STUDIES

In a descriptive, prospective, blinded image interpretation study⁶ of adult patients with known or suspected coronary artery disease, myocardial perfusion deficits (stress and rest PET images obtained with ammonia M 13 (n = 111) or rubidium Rb 82 chloride (n = 82) were compared. Changes in stenosis flow reserve (SFR) as determined by coronary angiography, PET perfusion defects at rest and stress for seven cardiac regions (anterior, apical, anteroposterior, posterolateral, posterolateral, and inferior walls) were graded on a scale of 0 (normal) to 5 (severe). Values for increasing flow reserve, defined as flow at maximum coronary vasodilatation relative to rest flow, ranged from 0 (total occlusion) to 5 (normal). With increasing impairment of flow reserve, the subjective PET defect severity increased. A PET defect score of 2 or higher was positively correlated with flow reserve impairment (SFR < 3).

A system review of published literature was conducted using pre-defined inclusion/exclusion criteria, which resulted in identification of 10 studies evaluating the use of Rb 82 PET myocardial perfusion imaging (MPI) for the identification of coronary artery disease as defined by catheter-based coronary angiography. In these studies, the patient was the site of analysis and 50% stenosis was the inclusion criterion. The original study on coronary artery disease (CAD) of these 10 studies, 9 studies of specificity (including 3 studies with 100% specificity). A total of 100% sensitivity and 100% specificity were included in a meta-analysis. The use of meta-analysis in establishing performance characteristics is limited, particularly by the possibility of publication bias (positive results being more likely to be published than negative results) which is difficult to detect especially when based on a limited number of small studies.

15 REFERENCES
 1. Senthimzichehian S. et al. Human biodistribution and radiation dosimetry of ⁸²Rb. J Nucl Med, 2010; 15:1592-99.
 2. Senthimzichehian S. et al. Radiation dosimetry of ⁸²Rb in humans under pharmacologic stress. J Nucl Med 2011; 52:485-91.
 3. Eckerman, K. F. et al. User's Guide to the DCAL System. ORNL/TM-2001-190, Oak Ridge National Laboratory, Oak Ridge, TN, August, 2006.
 4. Lederer, M and Shirley, V. Table of Isotopes, 7th Edition.
 5. Judge, S. et al. Applied radiation and isotopes (1987); vol 38, no. 3, pp 186-90.
 6. Demer, L.L. et al. Assessment of coronary artery disease severity by PET: Comparison with quantitative arteriography in 193 patients. Circulation 1989; 79:825-35.

16 HOW SUPPLIED/STORAGE AND HANDLING

16.1 How Supplied
 CardioGen-82[®] (rubidium Rb 82 generator), used to produce rubidium Rb 82 chloride injection, consists of Sr 82 adsorbed on a hydrous stannic oxide column with an activity of 3,330 MBq (90 millicuries) to 150 millicuries) Sr 82 at calibration time. A lead shield surrounded by a labeled plastic container encases the generator. The container label provides complete assay data for each generator. Use CardioGen-82 (rubidium Rb 82 generator) only with an appropriate, properly calibrated infusion system labeled for use with the generator.

16.2 Storage and Handling
 • Store the generator at 20° to 25°C (68° to 77°F) [See USP] Receipt, transfer, possession, storage, disposal or use of this product is subject to the radioactive material regulations and licensing requirements of the U.S. Nuclear Regulatory Commission (NRC), Agreement States or Licensing States as appropriate. Do not dispose of the generator in regular refuse systems.
 • For questions about the disposal of the CardioGen-82 generator, contact Bracco Diagnostics Inc. at 1-800-447-6883, option 3.

17 PATIENT COUNSELING INFORMATION

Pregnancy
 Advise a pregnant woman of the potential risk to a fetus *[see Use in Specific Populations (8.1)]*.

Lactation
 Advise lactating women that exposure to Rb-82 chloride through breast milk can be minimized if breastfeeding is discontinued when Rb 82 chloride injection is administered. Advise lactating women not to resume breastfeeding for at least one hour after completion of rubidium Rb 82 infusion. *[See Use in Specific Populations (8.2)].*

Post-Study Voiding
 Instruct patients to void after completion of each image acquisition session and as often as possible for one hour after completion of the PET scan. *[See Warnings and Precautions (5.5)].*

Manufactured for
 Bracco Diagnostics Inc.
 Monroe Twp, NJ 08831
 By GE Healthcare
 Merit-Physics, Inc.
 South Plainfield, NJ 07080

CardioGen-82[®] Rubidium Rb 82 Generator, Infusion Systems, and Accessories contain proprietary technology covered by one or more patents listed at www.braccodiagnosing.com/us-en/patents.

43-8200H Revised October 2020

g CardioGen·82[®]
(Rubidium Rb 82 Generator)

At the heart of cardiac PET

THE BASICS

Cardiogen-82 and the Model 1701 Infusion System

LIFE FROM INSIDE

ABBREVIATED PRESCRIBING INFORMATION

GL-CG-2300008 / April 2023

Contents

- 3** Introduction
- 4** What is Cardiogen-82?
- 5** What is the Model 1701 Infusion System?
- 6** Imaging with Cardiogen-82
- 7** Sample protocol for cardiac PET/CT imaging with rubidium-82 (⁸²Rb) chloride
- 8** Why choose cardiac PET with Cardiogen-82 and the Model 1701 Infusion System
- 12** Cardiogen Concierge Service
- 13** At the heart of cardiac PET
- 14** Bracco – your partner in cardiac PET
- 15** Abbreviated Prescribing Information

The use of cardiac PET has been growing during the last few years, with PET/CT scanners becoming widely available, and there is increasing evidence to demonstrate the advantages of this imaging approach. The publication of guidelines by the European Association of Nuclear Medicine (EANM) and, jointly, the Society of Nuclear Medicine and Molecular Imaging (SNMMI) and American Society of Nuclear Cardiology (ASNC) aims to standardise and support its use.^{1,2}

What is Cardiogen-82?

Cardiogen-82 is a radionuclide generator of rubidium-82 (⁸²Rb).

It is a closed system used to produce rubidium-82 (⁸²Rb) chloride injection for intravenous administration. Injection of rubidium-82 (⁸²Rb) chloride is indicated in PET imaging of the myocardium at rest and under pharmacological stress, in order to assess the regional myocardial perfusion in adults with a known or suspected coronary artery pathology.³


INDICATION

This medicinal product is for diagnostic use only.

The generator eluate (solution of rubidium-82 (⁸²Rb) chloride for injection) is used for Positron Emission Tomography (PET) imaging of the myocardium at rest or under pharmacologic stress conditions to evaluate regional myocardial perfusion in adults with suspected or known coronary artery disease.³

PET
Positron Emission Tomography

PET/CT
Positron Emission Tomography/Computed Tomography

**Generator
and cart
not shown
to scale**

What is the Model 1701 Infusion System?

The Cardiogen-82 Infusion System Model 1701 allows measurement and injection of precise activities of rubidium-82 (⁸²Rb) chloride that should not exceed 2200 MBq in unique activity, and 4400 MBq in cumulative activity, infused at a maximum rate of 50 mL/minute, with a maximum volume per infusion of 100 mL, and a total volume not exceeding 200 mL. Those limits for a session at rest and under stress reflect conditions for use in which clinical trials have been performed. Effective dose = 6 mSv when the maximum recommended cumulated radioactivity of 4440 MBq is administered.³

The Cardiogen-82 generator should only be used with the Cardiogen-82 Infusion System Model 1701, which has been designed specifically for use with the Cardiogen-82 generator.

Imaging with Cardiogen-82³

Patient information to be shared prior to scan

Pre-scan

- Fast for at least 6 hours prior to scan. Be well hydrated before the scan
- No caffeine for 12 hours prior to scan
- Women of child-bearing potential should have pregnancy excluded prior to the examination

Post-scan

- No particular precautions need to be taken following the administration of the radioactivity, but patients should be advised to urinate frequently in the first hour post-scan

Contraindications


- Pregnancy
- Children: safety and efficacy of rubidium-82 (⁸²Rb) chloride have not been established

Warnings

Strict adherence to instructions regarding quality control (QC) at clinical sites is required, with daily QC tests being performed to ensure patient safety.

The correct eluent for the generator must be used. Only a NaCl 9 mg/mL solution that meets pharmacopeia quality standards should be used for generator elution. The use of any other solutions is strictly forbidden as it may compromise patient safety.³

Sample protocol for cardiac PET/CT imaging with rubidium-82 (⁸²Rb) chloride³

APPROXIMATE TIMINGS


Note: In order to avoid the presence of residual activity from the previous dose, allow at least 10 minutes after completion of the previous Rb-82 infusion before starting the next dose.

Start imaging 70-90 seconds after the end of the injection of the Rb-82 injection in patients with a LVEF > 50%; if a longer circulation time is anticipated, the start of the acquisition may be delayed: 90-110 seconds if the LVEF is from 30% to 50%, 110 to 130 seconds if the LVEF is severely decreased, < 30% (i.e. a patient with severe malfunction of the left ventricle). In general, images acquisition lasts 5 minutes.^{2,3}

LVEF

Left Ventricular Ejection Fraction

**Generator
and cart
not shown
to scale**

Why choose cardiac PET with Cardiogen-82 and the Model 1701 Infusion System

**With its short half life of 75 seconds, you can provide
a high-quality patient-centred cardiac PET service with
Cardiogen-82^{1,2}**

- Fast rest/stress protocols (30 min) enabling high patient throughput^{1,3}
- Reduced radiation burden for both staff and patients compared to SPECT myocardial perfusion imaging (MPI)^{2,4}

Improve your image quality and diagnostic performance compared to SPECT MPI with Cardiogen-82¹

- High spatial and contrast resolution leads to high image quality⁵
- Cardiac PET MPI ⁸²Rb is superior to SPECT MPI in the detection of obstructive coronary artery disease (CAD)⁶
- Routine accurate, reliable attenuation correction decreases false positives, thus increases specificity^{2,5}
- Improved image quality and CAD detection in women and obese patients^{1,7}
- Ability to calculate myocardial blood flow (MBF) for both rest and stress allows assessment of coronary microvascular disease (CMD) which conveys additional diagnostic and prognostic information^{5,8}

The Model 1701 Infusion System is designed with you in mind

- Focused on efficiency, safety and ease of use
- Choice of flow rates
- Flex dosing without recalibration
- System self-monitoring with user alerts
- On-board record and report storage with re-print function, saving time and eliminating manual record keeping

Ready when you are – the Model 1701 Infusion System helps to make your day efficient

- On-board daily quality control with scheduling option to be run overnight, saving you time at the beginning of each day
- With no need for daily calibration, the system is ready to start the day with you
- Post-infusion flush makes the most of every dose and reduces the radiation exposure for users

With you all the way

- Support of the dedicated applications specialist team (the Cardiogen Concierge Service), provides you with the nuclear medicine expertise needed to start and maintain a successful cardiac PET programme

Cardiogen Concierge Service

With a dedicated team of accessible, responsive and knowledgeable expert applications specialists available to you in your time zone, the Cardiogen Concierge Service offers you support to start and maintain a successful cardiac PET service.


Cardiogen Concierge Service

The **Cardiogen Concierge Service** includes:

- Dedicated expert applications specialists with a named point of contact
- Support onsite in person for you and your entire nuclear medicine department
- Pre-training to introduce your team to Cardiogen-82 and the Model 1701 Infusion System
- Onsite ‘cold’ training with the Model 1701 Infusion System prior to first Cardiogen-82 generator delivery
- Presence for installation and first ‘hot’ patients
- Return for first generator exchange
- Certification of users with annual repeat training
- Performance of annual preventative maintenance checks
- Facilitation of inter-site relationships

With rapid response times, the Cardiogen Concierge Service is with you all the way, from service concept to implementation.

At the heart of cardiac PET

With the benefits of:

- Rapid throughput^{1,3}
- Superior diagnostic accuracy compared to SPECT MPI^{1,2,5,6}
- Low radiation exposure to both staff and patients^{2,4}
- Full support of the Cardiogen Concierge Service

...put the Cardiogen-82 generator and the Model 1701 Infusion System at the heart of your cardiac PET service

Bracco – your partner in cardiac PET

- Bracco has been at the forefront of cardiac PET for over 30 years and counting
- More than three million cardiac PET studies have been performed using Cardiogen-82 since 1989⁹
- Accessible, responsive and knowledgeable, Bracco professionals collectively bring more than 500 years of nuclear medicine expertise to help you start and maintain a successful cardiac PET service

Abbreviated Prescribing Information

References:

1. Sciegna R et al. *Eur J Nucl Med Mol Imaging*. 2021;48:1040-1069
2. Disizian V et al. *J Nucl Cardiol*. 2016;23:1187-1226
3. Cardiogen-82 (rubidium-82 generator (82Rb)) Common-EN SmPC (MRP) dated 25 November 2022
4. Schleipman AR et al. *J Nucl Cardiol* 2006; 13: 378-84
5. Schindler TH et al. *J Nucl Med*. 2020;Vol 6, 8:1221-1265
6. McArdle BA et al. *JACC*. Vol. 60, No. 18, 2012
7. Hyafil F et al. *J Nucl Cardiol*. 2020 Jun; 27(3): 755-768
8. Murthy VL et al. *J Nucl Cardiol*. 2018 Feb; 25(1): 269-297
9. Bracco sales data

CONFIDENTIAL

Cardiogen-82 Infusion System, Model 1701

Product	Notebook Number
Cardiogen-82	Not Applicable

Author(s):

Name	Function	Signature	Date
[REDACTED]	Principal Project Manager, Technical Operations	<i>Electronically signed-off</i>	01-Feb-2022

Reviewed by:

Name	Function	Signature	Date
[REDACTED]	Director, Pharmaceutics Operations	<i>Electronically signed-off</i>	01-Feb-2022
[REDACTED]	Associate Director, Technical Operations	<i>Electronically signed-off</i>	02-Feb-2022

Approval :

Name	Function	Signature	Date
[REDACTED]	Director, Device Regulatory Affairs	<i>Electronically signed-off</i>	10-Feb-2022
[REDACTED]	Director, Device Quality	<i>Electronically signed-off</i>	10-Feb-2022

Issuing date	14-Feb-2022
---------------------	-------------

Third party approval received.
 Name of Company (as applicable): Not Applicable

History of Changes

Previous Edition	Date	Reason(s) for Change
N/A		First Issuance

Table of Contents

History of Changes	2
1 Product Description	4
2 Intended Use.....	4
3 Usage Lifetime	4
4 Packaging	4
5 Sterilization	5
6 Quality Controls	5
7 Manufacturing and Conformity.....	5
8 Product Classification	5
9 Waste Disposal	5
10 Storage	5
11 Shelf Life.....	5
12 Age Stability.....	6

1 Product Description

The CardioGen-82[®] Infusion System Model 1701 is a reusable, electromechanical device contained within a mobile cart. Primary components of the system are the cart frame, lead radiation shielding, syringe pump, electronics including power supplies and printed circuit board assemblies (PCAs), and the enclosure. The lead radiation shielding is for shielding of the CG-82 Generator, syringe pump, positron detector, sterile disposable tubing set and components, and waste bottle. Radiation shielding also exists between the electronics and the rest of the primary cart components. The Na Constancy Source and Shield are mounted in the cart.

An on-board Panel PC with a graphical user interface (GUI), and software, provides integrated safety controls including Infusion System Calibration, and automated daily quality control procedures. The system software prevents patient infusions if specified quality control parameters or limits are not met, and system enforced limits for dose and volume is reached.

2 Intended Use

CardioGen-82[®] (rubidium Rb 82 generator) is a closed system used to produce rubidium Rb 82 chloride injection for intravenous administration. Rubidium Rb 82 chloride injection is indicated for Positron Emission Tomography (PET) imaging of the myocardium under rest or pharmacologic stress conditions to evaluate regional myocardial perfusion in adult patients with suspected or existing coronary artery disease.

The CardioGen-82[®] Infusion System is intended to be used with the CardioGen-82[®] Generator to deliver the radiopharmaceutical eluate rubidium (Rb 82) chloride injection to a patient in a controlled manner. Users of the CardioGen-82[®] Infusion System must be Nuclear Medicine Technologists. The CardioGen-82[®] Infusion System is intended to be used in a healthcare facility.

3 Usage Lifetime

The CardioGen-82[®] Infusion System Model 1701 enforces preventive maintenance on a yearly basis. As needed, parts are replaced or repaired to continue system operation.

4 Packaging

The Cardiogen-82[®] Infusion System, Model 1701 is secured inside a wooden crate for shipping. Trained service personnel uncrate and install the Cardiogen-82[®] Infusion System Model 1701 for use at each site.

5 Sterilization

The CardioGen-82[®] Infusion System Model 1701 is a non-sterile device.

6 Quality Controls

This device meets principle requirements and is in conformity with the relevant sections of the applicable harmonized EC standards and other normative documents (MDD 93/42/EEC and applicable recognized standards).

7 Manufacturing and Conformity

Bracco Diagnostic Inc. is the labeled manufacturer of the CardioGen-82[®] Infusion System, Model 1701.

8 Product Classification

The CardioGen-82[®] Infusion System Model 1701 is classified Class IIb accordingly to the MDD/93/42/EEC.

9 Waste Disposal

At the end of their useful lifetime, return CardioGen-82[®] Infusion System Model 1701 parts to authorized service personnel for recycling or disposal in accordance with local requirements.

10 Storage

The CardioGen-82[®] Infusion System should be stored at room temperature as noted on package labeling and in Operator Manual.

11 Shelf Life

The CardioGen-82[®] Infusion System Model 1701 has a ten (10) year defined shelf life. This product is not a single-use device.

12 Age Stability

The Cardiogen-82[®] Infusion System Model 1701 has no age stability criteria. Quality control is completed upon installation and every day before use.

CardioGen-82 NextGen2 System Exposure Readings
Old System vs New with Vulcan Lead, Day 0 and Day X Generator

Date: 12/15/15
 Generator Lot # 09101-15349X
 Day 0
 Potency on 12/15/15
 Sr-82 Activity: 125 mCi
 Sr-85 Activity: 95 mCi

Date: _____
 Generator Lot # _____
 Day _____
 Potency on _____
 Sr-82 Activity: _____ mCi
 Sr-85 Activity: _____


All measurements in mR/hr

Location (Reference attached cart location orientation)	Day 0 Original Cart	Day 0 New Cart	Day __ Original Cart	Day __ New Cart
Front- shield surface	6.500	6.500		
Front- cart body surface	2.500	3.000		
Front- 1 meter from cart surface	.070	.070		
Left- shield surface	Not accessible	6.500	Not accessible	
Left- cart body surface	1.000	2.500		
Left- 1 meter from cart body surface	.060	.060		
Right- shield surface	Not accessible	6.200	Not accessible	
Right- cart body surface	1.300	2.500		
Right- 1 meter from cart body surface	.060	.060		
Back- cart body surface	.020	.300		
Back- 1 meter from cart body surface	.006	.025		
Top- lid surface	1.500	3.500		
Top- cart body	.350	.500		
Bottom surface	14.000	1.500		
1 foot halo from side of cart body (max reading, location)	.400	.400		
1 meter halo from side of cart body (max reading, location)	.070	.070		
Background	.006	.006		
Control Panel (approx. waist high)	.090	.090		


Generator Day 0 Surface = 180 mR/hr, 1 foot = 9.0 mR/hr, 1 meter = 1.0 mR/hr
 Generator Day __ Surface = _____ mR/hr, 1 foot = _____ mR/hr, 1 meter = _____ mR/hr

Survey Instrumentation

Bicron, micro rem, RbM MTI# 0712 cal date: 3/32/15
 Note- these are checked daily with a calibrated source.

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Product identifier Sheet Code: 279

Trade name: CardioGen-82

Chemical Name: Rubidium Rb-82 generator for elution of rubidium chloride Rb-82 injection.

Synonyms: For active: Rubidium Chloride 82, Rubidium Chloride Rb-82 Injection.

CAS Number: 7791-11-9

Application of the substance / the mixture:

A myocardial perfusion agent useful in distinguishing normal from abnormal myocardium.

We recommend that you use this product in a manner consistent with the listed use. If your intended use is not consistent with the stated use, please contact your sales or technical service representative.

Rubidium Rb-82 generator supplied in the form of strontium Sr-82 adsorbed on a hydrous stannic oxide column which is encased in a lead shield and surrounded by a plastic container.

Details of the supplier of the safety data sheet

Manufacturer/Supplier:

Bracco Diagnostics Inc.

P.O. Box 5225

Princeton, NJ

08543

Phone number: 1-800-257-5181

Email: HSE@bracco.com (responsible for the SDS)

Information department:

B-Lands Consulting

WTC, 5 Place Robert Schuman, BP

1516 38025 Grenoble, FRANCE

Tel: +33 476 295 869

Fax: +33 476 295 870

Email: clients@reachteam.eu

www.reachteam.eu

Emergency telephone number:

EMERGENCY CONTACT:

Health: 1-800-257-5181

U.S. Transport - Chemtrec: 1-800-424-9300

International Transport - Chemtrec: 1-703-527-3887

2 Hazard(s) identification

Classification of the substance or mixture

The product is not classified, according to the Globally Harmonized System (GHS).

Label elements

GHS label elements Not applicable.

Hazard pictograms Not applicable.


Signal word Not applicable.

Hazard statements Not applicable.

Additional Information:

Classification system:

NFPA ratings (scale 0 - 4)

Health = 0

Fire = 0

Reactivity = 0

Safety Data Sheet

According to OSHA HazCom Standard [2012]

Printing date 12/12/2018

Version 4

Reviewed on 03/05/2022

HMIS-ratings (scale 0 - 4)

HEALTH	0	Health = 0
FIRE	0	Fire = 0
REACTIVITY	0	Reactivity = 0

Other hazards No further relevant information available.

Results of PBT and vPvB assessment

PBT: Not applicable.

vPvB: Not applicable.

3 Composition/information on ingredients

Chemical characterization: Mixtures

Description: Mixture: consisting of the following components.

Hazardous Components			
CAS No.	Name	Classification	Qty.
1310-73-2	Sodium Hydroxide	Skin Corr. 1A, H314, Acute Tox. 4, H302	2 – 2.5 %
7440-24-6	Strontium ⁸² / Strontium ⁸⁵ (**)	Water-react. 1, H260	0.1 – 1 %

Non-Hazardous components		
CAS No.	Name	Qty.
18282-10-5	Tin (IV) Oxide	1 – 10 %

Information on components		
CAS No.	Name	Qty.
18282-10-5	Tin (IV) Oxide	1 – 10 %

Additional information:

** CAS number pertains to Strontium

When sterile, pyrogen-free Sodium Chloride Injection USP is used to elute the generator, the diagnostic agent Rubidium Chloride Rb-82 injection is formed. The resulting material is radioactive.

Rubidium Rb-82 (eluent from CardioGen-82 generator) decays by positron emission and associated gamma emission with a physical half-life of 75 seconds. The specific gamma ray constant for Rb-82 is 6.1 R /h-mCi at 1 cm. The first half-value layer is 0.7 cm of lead (Pb).

4 First-aid measures

Description of first aid measures

General information:

Formulation contains RADIOACTIVE CHROMIUM 51.

As in the use of any radioactive material, care should be taken to insure minimum radiation exposure. Personnel who handle radioactive materials should be trained in their use and should follow appropriate precautions for work with these materials.

See Health Effects and Toxicology sections for additional information.

After inhalation:

Supply fresh air. If required, provide artificial respiration. Consult doctor if symptoms persist.

After skin contact:

Immediately wash with water and soap and rinse thoroughly. Seek medical attention if skin irritation, swelling or redness develops and persists. Swab skin to verify all removable radioactive contamination is gone.

Safety Data Sheet

According to OSHA HazCom Standard [2012]

Printing date 12/12/2018

Version 4

Reviewed on 03/05/2022

After eye contact:

Rinse immediately with plenty of water, also under the eyelids for at least 15 minutes. If irritation persists get medical attention. Swab eyelid and surrounding area to verify all removable radioactive contamination is gone.

After swallowing:

Immediately call a doctor.

Vomiting may be induced only if a person is conscious and if ingestion has occurred within the past three hours. Never induce vomiting in a person who is unconscious or experiencing convulsions. If ingestion of eluent from the CardioGen-82 generator (containing radioactive Rubidium-82) inadvertently occurs, the individual may be treated by water hydration or diuresis to facilitate elimination of the radioactive material.

Most important symptoms and effects, both acute and delayed

No further relevant information available.

Indication of any immediate medical attention and special treatment needed

No further relevant information available

5 Fire-fighting measures

Extinguishing media

Suitable extinguishing agents: Water

Special hazards arising from the substance or mixture

WARNING: Material is RADIOACTIVE

Advice for firefighters

Evacuate personnel to an upwind direction, remove unneeded material and cool container(s) with water from a maximum distance. Move container from fire area if you can do it without risk.

Protective equipment:

Firefighters should wear adequate personal protective equipment with protection of respiratory tract (self-contained breathing apparatus) (SCBA). In addition, firefighters should wear flame and chemicals resistant clothing, boots and gloves.

6 Accidental release measures

Personal precautions, protective equipment and emergency procedures

Wear protective equipment. Keep unprotected persons away.

Radioactive shielding if appropriate to minimize exposure

Environmental precautions:

Do not allow product to reach sewage system or any water course. Do not allow to penetrate the ground/soil.

Inform respective authorities in case of seepage into water course or sewage system.

Methods and material for containment and cleaning up:

Absorb with liquid-binding material.

Place spilt material in an appropriate container for disposal.

The spill area should be ventilated and decontaminated after material is collected.

Dispose of radioactive material in accordance with all local, state, federal and EC Regulations or with the Regulations of the country in which the material is used.

Reference to other sections

No dangerous substances are released.

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

Protective Action Criteria for Chemicals

PAC-1:

None of the ingredients are listed.

Safety Data Sheet

According to OSHA HazCom Standard [2012]

Printing date 12/12/2018

Version 4

Reviewed on 03/05/2022

PAC-2:

None of the ingredients are listed.

PAC-3:

None of the ingredients are listed.

7 Handling and storage

Precautions for safe handling

Avoid contact with the eyes and skin.

Radioactive materials may pose significant health risks if not properly handled. Personnel who handle radioactive materials should be trained in their use and should follow appropriate precautions for work with these materials.

Information about protection against explosions and fires:

No special measures required.

Conditions for safe storage, including any incompatibilities

Requirements to be met by storerooms and receptacles:

Store in a cool, dry place in well-sealed receptacles.

Container Requirements: The CardioGen-82 (Rubidium Rb 82 Generator) is encased in a lead shield surrounded by a labeled plastic container.

Storage Conditions: Store at 20-25 degrees C.

Information about storage in one common storage facility:

Not required.

Further information about storage conditions:

None.

Specific end use(s)

No further relevant information available.

8 Exposure controls/personal protection

Control parameters

Components with limit values that require monitoring at the workplace:

The following constituent is the only constituent of the product which has a PEL, TLV or other recommended exposure limit.

At this time, the remaining constituent has no known exposure limits.

1310-73-2 Sodium Hydroxide

TLV-ACGIH	Long-term value: 2 mg/m ³ absolute limit value for 2009
-----------	-----------------------------------------------------------------------

Additional information:

The lists that were valid during the creation were used as basis.

The Exposure Limit for Stannic Oxide as Tin, is 2 mg / m³ for the OSHA Permissible Exposure Limit (PEL) and for the American Conference of Governmental Hygienists (ACGIH) Threshold Limit Value-Time Weighted Average (TLV-TWA).

The Occupational Guideline to Exposure to Radioactive Materials is an effective dose equivalent to 50 mSV per year.

Exposure controls

Appropriate Technical Controls:

Provide adequate aspiration / ventilation in the workplace

Radioactive shielding if appropriate to minimize exposure

Personal protective equipment

General protective and hygienic measures:

The usual precautionary measures for handling chemicals should be followed. Wash hands before breaks and at the end of work.

Safety Data Sheet

According to OSHA HazCom Standard [2012]

Printing date 12/12/2018

Version 4

Reviewed on 03/05/2022

Do not eat, drink and smoke while working.

Breathing equipment:

Not anticipated for normal clinical environment.

In non-routine exposure conditions, where risk assessment shows air-purifying respirators are appropriate, use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Self-contained breathing apparatus should be available for emergency use.

Protection of hands:

Protective gloves

The glove material must be impermeable and resistant to the product/ the substance/ the preparation. Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

Material of gloves

Natural rubber, NR Nitrile rubber, NBR

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material cannot be calculated in advance and has therefore to be checked prior to the application.

Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

Eye protection:

Safety glasses

Body protection: Protective work clothing.

9 Physical and chemical properties

Information on basic physical and chemical properties

General Information:

Appearance

Form:	Solid
Color:	Not applicable
Odor:	Not determined
Odor threshold:	Not determined
pH – Value:	Not applicable
Melting point/Melting range:	Not determined
Boiling point/Boiling range:	Not determined
Flash point:	Not applicable.
Flammability (Solid, Gaseous):	Not determined
Ignition temperature:	Product is not self-igniting

Safety Data Sheet

According to OSHA HazCom Standard [2012]

Printing date 12/12/2018

Version 4

Reviewed on 03/05/2022

Danger of explosion:	Product does not present an explosion hazard
Flammability Limits	
Lower:	Not determined
Upper:	Not determined
Explosion Limits	
Lower:	Not determined
Upper:	Not determined
Oxidizing properties:	Not determined
Vapor pressure:	Not applicable
Density:	Not determined
Relativity density at 20 °C (68 °F):	Not determined
Vapor Density at 20 °C (68 °F):	Not determined
Evaporation rate:	Not determined
Solubility in / Miscibility with	Not determined
Water:	Insoluble
Partition coefficient (n-octanol/water):	Not determined
Viscosity	
Dynamic:	Not determined
Kinematic:	Not determined
Water:	Not determined
VOC Content:	Not determined
Solid content:	11.0 %

Other Information

No further information available

10 Stability and reactivity

Reactivity

There are not particular dangerous reactions with other substances in normal conditions of use.

Chemical stability

This generator contains strontium Sr-82 and Sr-85. Sr-82 decays with a half-life of 25 days and Sr-85 with a half-life of 65 days.

Rubidium Rb-82 (eluent from CardioGen-82 generator) decays by positron emission and associated gamma emission with a physical half-life of 75 seconds. The specific gamma ray constant for Rb-82 is 6.1 R/h-mCi at 1 cm. The first half-value layer is 0.7 cm of lead (Pb).

Thermal decomposition / conditions to be avoided:

No decomposition if used according to specifications.

Possibility of hazardous reactions No dangerous reactions known.

Conditions to avoid No further relevant information available.

Incompatible materials: No further relevant information available.

Hazardous decomposition products: No dangerous decomposition products known.

Safety Data Sheet

According to OSHA HazCom Standard [2012]

Printing date 12/12/2018

Version 4

Reviewed on 03/05/2022

11 Toxicological information

Information on toxicological effects

Acute toxicity:

LD/LC50 values that are relevant for classification:		
18282-10-5 TIN (IV) OXIDE		
Oral	LD50	> 20000 mg/kg (Rat)
7647-14-5 sodium chloride		
Oral	LD50	3000 mg/kg (Rat)
1310-73-2 Sodium Hydroxide		
Oral	LD50	2000 mg/kg (Rat)
hydrochloric acid		
Oral	LD50	900 mg/kg (Rabbit)

Primary irritant effect:

on the skin:

No irritant effect.

Material contains low concentration of components that are mild irritants or possible irritants. It may have potential to cause mild irritation, however, moderate or severe irritation is not expected. Acute skin exposure is not expected to exceed 50 mSV/year (Occupational Exposure Guideline)

on the eye: No irritating effect.

Sensitization:

No sensitizing effects known. Biological effects from exposure to radioactive materials are based on exposures higher than those permitted in an occupational setting. No harmful effects are expected from CardioGen-82 under normal use conditions

Other information (about experimental toxicology):

By Inhalation:

Formulation contains some materials that are irritants. Inhaling small amounts of liquid aerosol may result in nasal and other respiratory tract irritation. In general, depending on both the dose and duration, exposure to radioactive materials may produce adverse effects.

By Ingestion:

Inadvertent ingestion of trace amounts of this material would not be expected to result in symptoms. However, ingestion of substantial amounts may cause kidney damage (acute tubular necrosis). In general, depending on both the dose and duration, exposure to radioactive materials may produce adverse effects. Toxic effects of radiation may include development of cancer, adverse reproductive effects, including adverse effects on the fetus.

Additional toxicological information:

The product is not subject to classification according to internally approved calculation methods for preparations. When used and handled according to specifications, the product does not have any harmful effects to our experience and the information provided to us. Contact with small quantities of material for short periods is not expected to result in pharmacologic or toxic effects.

Carcinogenic categories

IARC (International Agency for Research on Cancer)
None of the ingredients are listed.
NTP (National Toxicology Program)
None of the ingredients are listed.
OSHA-Ca (Occupational Safety & Health Administration)
None of the ingredients are listed.

Safety Data Sheet

According to OSHA HazCom Standard [2012]

Printing date 12/12/2018

Version 4

Reviewed on 03/05/2022

12 Ecological information

Toxicity

Aquatic toxicity:	
18282-10-5 TIN (IV) OXIDE	
EC50/48 h	> 100 mg/l (Daphnia Magna)
NOEC/48h	> 100 mg/l (Daphnia Magna)
1310-73-2 Sodium Hydroxide	
LC50	180 mg/l (Fish)

Persistence and degradability

No further relevant information available.

Bio accumulative potential

No further relevant information available.

Mobility in soil

No further relevant information available.

Results of PBT and vPvB assessment

PBT: Not applicable.

vPvB: Not applicable.

Other adverse effects

No further relevant information available.

13 Disposal considerations

Waste treatment methods

Recommendation:

Must not be disposed of together with household garbage. Do not allow product to reach sewage system. Reutilize if possible or contact a waste processor for recycling or safe disposal.

Uncleaned packaging:

Recommendation:

Dispose of radioactive material in accordance with all local, state, federal and NRC regulations or with the regulations of the country in which the material is used.

14 Transport information

UN – Number

DOT:

1. Domestic

DOT Shipping Classification: Radioactive Material.

See US Regulations regarding shipment of radioactive materials.

2. International

DOT Shipping Classification: Radioactive material.

See US and other Nations Regulations on shipment of radioactive materials.

ADR:

1. Domestic

ADR Shipping Classification: Radioactive Material.

See US Regulations regarding shipment of radioactive materials.

2. International

ADR Shipping Classification: Radioactive material.

See US and other Nations Regulations on shipment of radioactive materials.

Safety Data Sheet

According to OSHA HazCom Standard [2012]

Printing date 12/12/2018

Version 4

Reviewed on 03/05/2022

IMDG:

1. Domestic

IMDG Shipping Classification: Radioactive Material.

See US Regulations regarding shipment of radioactive materials.

2. International

IMDG Shipping Classification: Radioactive material.

See US and other Nations Regulations on shipment of radioactive materials.

IATA:

1. Domestic

ICAO/IATA Shipping Classification: Radioactive Material.

See US Regulations regarding shipment of radioactive materials.

2. International

ICAO/IATA Shipping Classification: Radioactive material.

See US and other Nations Regulations on shipment of radioactive materials.

UN proper shipping name

DOT Not applicable.

Transport hazard class(es)DOT

Class

Not applicable.

Label

8.0

Packing group

DOT

Not applicable.

Environmental hazards:

Marine pollutant:

No

Special precautions for user:

Not applicable.

Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code

Not applicable.

**15 Regulatory information****Safety, health and environmental regulations/legislation specific for the substance or mixture**

Sara

Section 355 (extremely hazardous substances):	
None of the ingredients are listed.	
Section 313 (Specific toxic chemical listings):	
None of the ingredients is listed.	
TSCA new (21st Century Act) (Substances not listed)	
1310-73-2	Sodium Hydroxide
7440-24-6	Strontium ⁸² /Strontium ⁸⁵ (**)

Proposition 65

Chemicals known to cause cancer:
None of the ingredients is listed.
Chemicals known to cause reproductive toxicity for females:
None of the ingredients are listed.
Chemicals known to cause reproductive toxicity for males:
None of the ingredients are listed.

Safety Data Sheet

According to OSHA HazCom Standard [2012]

Printing date 12/12/2018

Version 4

Reviewed on 03/05/2022

Chemicals known to cause developmental toxicity:

None of the ingredients are listed.

Carcinogenic categories

EPA (Environmental Protection Agency)

None of the ingredients are listed.

TLV (Threshold Limit Value established by ACGIH)

None of the ingredients are listed.

NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients are listed.

GHS label elements Not applicable.

Hazard pictograms Not applicable.

Signal word Not applicable.

Hazard statements Not applicable.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

Description and Property:

The CardioGen-82 Rubidium Rb 82 Generator is a shielded column that contains a slurry of normal saline solution and stannic oxide to which Strontium Sr-82/Sr-85 has been adsorbed.

As in the use of any radioactive material, care should be taken to insure minimum radiation exposure.

Significant Dangers:

If ingestion of eluent from the CardioGen-82 generator (containing radioactive Rubidium-82) inadvertently occurs, the individual may be treated by water hydration or diuresis to facilitate elimination of the radioactive material.

Training Hints:

All persons handling this product should be informed on the existence of the hazard, on any possible risk they might be subjected to and about all required protective measures to prevent such a damage or to reduce the exposition.

WARNINGS:

When transporting an employee for medical assistance, after the employee has had direct contact with a radioactive material, care should be taken to avoid contamination of transport vehicle and medical facility. Skin decontamination and monitoring should be conducted as appropriate.

Radiopharmaceuticals should be used only by physicians who are qualified by training and experience in the safe use and handling of radionuclides) and whose experience and training have been approved by the appropriate government agency authorized to license the use of radionuclides.

Dispose of radioactive material in accordance with all local, state, federal and NRC regulations or with the regulations of the country in which the material is used.

Contact:

Bracco Diagnostics Inc.

P.O. Box 5225

Princeton, NJ

08543

Date of preparation / last revision 03/05/2022, revision 4

Changes: General revision of the entire Safety Data Sheets, modified section 3, section 9, section 14 and section 16.

Abbreviations and acronyms:

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

Safety Data Sheet

According to OSHA HazCom Standard [2012]

Printing date 12/12/2018

Version 4

Reviewed on 03/05/2022

ACGIH: American Conference of Governmental Industrial Hygienists
EINECS: European Inventory of Existing Commercial Chemical
Substances ELINCS: European List of Notified Chemical Substances
CAS: Chemical Abstracts Service (division of the American Chemical Society)
NFPA: National Fire Protection Association (USA)
HMIS: Hazardous Materials Identification System (USA)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

Water-react. 1: Substances and mixtures which in contact with water emit flammable gases – Category 1

Acute Tox. 4: Acute toxicity – Category 4

Skin Corr. 1A: Skin corrosion/irritation – Category 1A

*** Data compared to the previous version altered.**

- data updating on the basis of the latest amendments.